



## MShorts 3 – Technical Specification

<b>Name</b>	MShorts 3 (part of Mbody 3 product family)
<b>Type</b>	#52-1882 MB ... #52-2282 MB
<b>Used with</b>	MCell 3
<b>For applications</b>	Mbody Live 3 (iOS 9 or later, Android 6 or later) Muscle Monitor 3
<b>Channels</b>	6 EMG biosignal channels (3 on left side, 3 on right side):  Quadriceps group Hamstrings group Gluteus
<b>Textile materials</b>	Carvico Revolutionale® 71% Polyamide micro, 29% Elastane (lycra xtralife) (195g/m2) Electrodes silver coated yarns
<b>Washing</b>	Fully machine washable at 30°C (86°F)   Avoid exposing sensors to solvent.
<b>Operating temperature</b>	0°C ... +40°C (32°F ... + 104°F)
<b>Humidity</b>	10 ... 90 % RH (non-condensing)
<b>Water resistance</b>	Protected against water splashes during use. Waterproof during washing. Let dry carefully before use with MCell 3.
<b>Sizes</b>	S, M, L, XL, XXL (Unisex)
<b>Weight</b>	~ 190 g (6,7 oz)
<b>Recycling</b>	 This product contains electronic components and materials. Product must be disposed as electrical and electronic equipment.
<b>Standards and approvals</b>	Oeko-Tex® Standard 100 CE, FCC, CMIIT, RoHS, WEEE